

A guide to using the new Stallergenes Greer



We're excited to offer this convenient tool to help you calculate and prepare dilutions of allergenic extracts for use in testing and treating your patients. The tool can be used for PNU/mL, weight/volume, AU/mL, and BAU/mL calculations.

Just go to [STAGReextractcalculator.com](https://stallergenesgreer.com) and enter your PIN on the homepage, and the main screen will appear



Another feature you'll appreciate:

- Ability to create a shortcut by saving the domain as an icon on your phone so you can easily access the calculator



To see for yourself, sign in at [STAGReextractcalculator.com](https://stallergenesgreer.com)



Please contact Medical Affairs at 855-752-5046 for additional questions regarding PNU/mL dilution calculations for glycerinated solutions using this calculator.

For questions regarding any of our other products or services, please call 800-378-3906 or email uscustomer@stallergenesgreer.com.

We're here for you!

DISCLAIMER: This website, including the dilution calculator, is for informational purposes only and is intended to be used by qualified healthcare professionals. All calculations should be confirmed before use. Any reliance placed on the content of the website is at your own risk. Stallergenes Greer makes no claims as to the accuracy of the information contained herein, and the information provided is not a substitute for the clinical judgment of a qualified healthcare professional. The calculator is not intended to be used for medical diagnosis or treatment, or as a substitute for a medical diagnosis and/or treatment rendered or prescribed by a qualified healthcare professional. Stallergenes Greer does not directly or indirectly practice medicine or provide medical services, and therefore assumes no liability whatsoever of any kind for the information produced by the calculator or for any diagnosis or treatment made in reliance thereon. Neither Stallergenes Greer nor any other party involved in the preparation or publication of this tool shall be liable for any decisions made or actions taken by you or others in reliance on this information.

Stallergenes Greer will not be collecting any personally identifiable information (PII). Please see additional information regarding our privacy practices at www.stallergenesgreer.com/privacy-policy/.

STALLERGENES  GREER

©Stallergenes Greer 2023. All Rights Reserved. 082823H2873 09/23


Beyond supplying quality extracts to your practice

We're putting a convenient new tool in your hands


Introducing the Stallergenes Greer



**Dilution
Calculator**
Dilute extracts in 3 easy steps




STALLERGENES  GREER




STEP 1

Enter the concentration of the bulk extract you wish to dilute by selecting either PNU/mL, weight/volume, AU/mL, or BAU/mL.



STEP 2

Enter your desired final strength by either selecting a preset value or entering a custom strength.



STEP 3

Select your desired final volume or enter a custom amount. Tap "Calculate."



Your final formula will appear, with mixing instructions.

After you've reviewed your final calculations, you can print them, and then clear them to make another calculation.